

High-performance quantitative exposomics covering up to >230 toxicants and key biomarkers

Yunyun Gu^{a,b}, Max L. Feuerstein^{a,c}, Caroline H. Johnson^d, Benedikt Warth^{a,b,c*}

^aDepartment of Food Chemistry and Toxicology, Faculty of Chemistry, University of Vienna, Vienna, Austria

^bVienna Doctoral School of Chemistry, University of Vienna, Währinger Straße 42, 1090 Vienna, Austria

^cExposome Austria, Research Infrastructure and National EIRENE Node, Austria

^dDepartment of Environmental Health Sciences, Yale School of Public Health, New Haven, CT 06510, United States of America

CORRESPONDING AUTHOR: Benedikt Warth, E-mail: benedikt.warth@univie.ac.at, Phone: +43 1 4277 70806, Department of Food Chemistry and Toxicology, Faculty of Chemistry, University of Vienna, Währinger Straße 38, 1090 Vienna, Austria.

KEYWORDS. *Next-generation human biomonitoring, Mass spectrometry, Early-life chemical exposure, Exposome wide association study (ExWAS), Public and environmental health*

ABSTRACT: The exposome encompasses environmental exposures throughout life and significantly impacts health and disease. Exposure chemicals, typically present at trace levels, are mostly quantified using targeted LC-MS/MS. However, many existing methods are limited to a narrow range of analyte classes or lack sufficient sensitivity for exposomic analyses and the applicability to large sample cohorts for exposome-wide association studies (ExWAS) remains to be demonstrated. Here, we present a scalable workflow for analyzing >230 biomarkers in urine, plasma, and serum using solid-phase extraction in 96-well plates and LC-MS/MS. Moreover, a new conceptual framework for validation criteria of assays designed to analyze highly diverse compounds at trace levels is proposed. Method robustness was evaluated, demonstrating suitable extraction recovery and matrix effects (SSE) within 60-130%, inter-/intra-day precision (RSD) <30%, and exceptional sensitivity (limit of detection, 0.015-50 pg/mL) for 60-80% of the analytes across the investigated matrices. To showcase the method's applicability in epidemiological studies, 200 urine samples from pregnant women in a longitudinal cohort were analyzed, with more than 130 biomarkers detected in the real-life samples, several for the first time in US urine. With its broad analyte coverage, excellent sensitivity, robustness, and exceptionally high sample throughput, this method offers the necessary performance for large-scale ExWAS studies in the future.

Individuals are exposed to thousands of chemicals through dietary intake and other environmental sources. Exposomics aims to comprehensively assess these environmental exposures and to identify potential links to adverse health effects, including chronic diseases or cancer (Koual et al., 2020; Wild, 2005). Liquid chromatography coupled to mass spectrometry (LC-MS) is currently the most relevant technique in exposome-wide association studies (ExWAS) and can be used for targeted and non-targeted analysis (Lai et al., 2024). Targeted assays focus on quantifying a pre-defined set of compounds with high sensitivity, accuracy, and precision (Jagani et al., 2022; Jamnik et al., 2022). In contrast, non-targeted analysis and suspect screening, aim to identify and characterize chemicals without a priori knowledge of the sample composition and can lead to discovery of novel or unexpected exposures (Manz et al., 2023). Due to favorable assay robustness, sensitivity and quantitative capabilities, multi-analyte targeted LC-MS/MS methods are often the preferred choice. This allows detection and quantification of a wide range of compounds at trace concentration levels and provides the basis for exposomics (Braun et al., 2020; Krausová et al., 2024). Suitable sample preparation protocols are essential, as sample

matrix components can substantially impact signal intensities, particularly for trace level compounds (Preindl et al., 2019). In addition, limited duty cycle can hamper data quality in large-scale LC-MS/MS assays targeting hundreds of analytes in a single run (Cai & Yan, 2021). Various sample preparation protocols, including protein precipitation (PPT) and solid phase extraction (SPE), are commonly used to extract exposure biomarkers from diverse human matrices and aim to provide best assay sensitivity and robustness, while reducing unwanted instrumental contaminations at the same time. PPT protocols are widely used in omics studies due to their broad analyte coverage (Huang et al., 2023; Jamnik et al., 2022), whereas SPE methods are traditionally tailored to extract and analyze a smaller number of analytes or analyte classes. Yet, SPE is increasingly popular for exposome-scale studies due to its efficient cleanup of complex biological samples (Gu et al., 2023). Recently, a generic SPE protocol for extracting 94 exposure compounds with diverse physicochemical properties was developed and demonstrated its applicability to targeted and non-targeted exposomics (Gu et al., 2025). Analyte coverage was compared to a generic PPT protocol and results showed a substantial overlap of

identified analytes of the PPT and SPE method in urine and plasma. Broad analyte coverage is critical for targeted multi-class analysis, but SPE protocols for >200 compounds remain rare. One of the key benefits of SPE methods is the capability to perform high-throughput sample extraction in 96-well plates. Time requirements and consumption of single-use labware substantially impacts the economic and environmental footprint of sample preparation, which is especially relevant for ExWAS where hundreds to thousands of samples need to be analyzed (Gu et al., 2023; Lin et al., 2022).

To ensure suitable data quality, rigorous method validation is essential. This is especially important when analytes are present at trace levels. The latest update of the guideline from European Commission (EC) No. 2021/808 provides systematic criteria (referred to as “EC criteria”) for method validation in terms of trueness (evaluated via extraction recovery, RE), precision (intermediate precision, RSDR), and repeatability (RSDr), over a wide range of concentrations (EC, 2021). However, the EC criteria only provide limited applicability for exposome-scale LC-MS/MS methods, because of the complex nature of large-scale multianalyte methods and the extremely low concentrations of most exposure compounds in human matrices. In particular, the EC criteria do not provide clear limits for precision and repeatability for analyte concentrations below 120 µg/kg and do not account for the large number of analytes in extreme multi-analyte methods required for exposomics. Therefore, new systematic validation criteria for large-scale exposomics assays for trace analysis are required to address these challenges and to further complement EC criteria. Finally, applicability for large-scale studies needs to be tested to showcase the robustness of exposomics workflows, which is of special relevance when for populations with low-level exposures like pregnant women are targeted (Krausová et al., 2024).

In this study, a multi-analyte LC-MS/MS method for 234 highly diverse chemicals covering a wide range of polarities (LogP between -4.6-9.6) with detection limits (LOD) in the pg/mL-range was developed. To enhance sample cleanup, assay sensitivity and sample throughput, a recently developed SPE protocol for exposomics in 96-well plates was used and the workflow was validated for common human sample matrices like urine, plasma, and serum. Based on empirical data extracted from relevant literature and generated within this study, we propose a systematic conceptual framework for method validation in exposomics. The new approach was used together with criteria set by the European Commission to assess fitness-for-purpose. Finally, validation results were used to systematically classify assay compounds in order to report results as quantitative, semi-quantitative, or qualitative. To assess the feasibility of the assay, it was applied to 200 urine samples obtained from pregnant US females during gestation (12th, 20th, 28th, and 36th week).

EXPERIMENTAL SECTION

Chemicals and materials. The diverse analyte panel consisted of 234 compounds and was developed based on previous work (Braun et al., 2022; Hossain et al., 2024; Jamnik et al., 2022; Preindl et al., 2019). Furthermore, 58 additional endocrine disrupting chemicals (EDCs) were selected based on relevant databases, i.e. compounds listed by the U.S. Environmental Protection Agency (EPA), European Human Biomonitoring Initiative (HBM4EU), Comparative Toxicogenomics Database (CTD), and relevant human biomonitoring (HBM) studies (see Table S1). In total, this panel included 13 classes of compounds including perfluorinated alkylated substances (PFAS), medical drugs, personal care products, air pollutants, pesticides, flame retardants, mycotoxins, industrial products, food processing by-products, plastics-related chemicals, disinfection and by-products, phytotoxins/-estrogens, and endogenous estrogens. Working stock solutions of all analytes (referred to as “STD mix”) were created by diluting individual stock solutions with acetonitrile. The following isotope-labeled

chemicals were used as internal standards: 13C12-bisphenol A, 13C6-butylparaben, 13C15-deoxynivalenol, 2H3-erythromycin, 13C6-ethylparaben, 13C2-mono-butyl phthalate, 13C6-methylparaben, 13C8-perfluorooctanoic acid, 13C8-perfluorooctanesulfonic acid, 13C6-propylparaben and 13C18-zearalenone. To account for variations in the SPE process and instrumental analysis, 13C18-zearalenone (100 ng/mL) was added before sample preparation, while all other internal standards (referred to as “ISTD mix”) were added after SPE to track instrumental performance and matrix effects. Used concentrations are provided in Table S2 and product numbers and suppliers are listed in Table S3. Oasis PRiME HLB 96-well plates (30 mg, 2 mL) were purchased from Waters. Further details regarding the preparation of stocks, mixtures, solvents, and the spiking process are presented in the Supplementary Information (SI).

Samples. Urine samples (n = 200) were collected from 50 pregnant women in Connecticut (U.S.) participating in the Yale Pregnancy Outcome Prediction Study (YPOPS). Participants self-collected samples at gestational weeks (GW) 12, 20, 28, and 36 during a single-day visit. Urine was collected into sterile cups, kept on ice, and sent to a collection center, where 1.5 mL aliquots were transferred to cryovials, stored at -80°C, and shipped on dry ice to the Global Exposomics and Biomonitoring Laboratory at the University of Vienna for analysis. Previously, these samples have been analyzed for biomarkers of mycotoxin exposure (Krausová et al., 2024). The study was approved by the Yale University Human Investigation Committee (#1601017004), and samples were banked by the Yale University Reproductive Sciences (YURS) Biobank (HIC#1309012696).

Pooled urine, plasma, and serum samples were used for method development and validation. The pooled urine was collected from a female volunteer following a two-day period of abstaining from foodstuff and beverages stored in plastic containers, phytoestrogen-rich foods, and cosmetics containing parabens (Preindl et al., 2019). Pooled plasma was purchased from Innovative Research (IPLAWBLIH-31982) (Novi, MI, USA), and serum from Sigma-Aldrich (H4522, Sigma-Aldrich, Germany). All samples were aliquoted and stored at -80 °C until analysis.

Sample preparation. An optimized SPE protocol utilizing 96-well plates was used according to a recently developed protocol (Gu et al., 2025). In short, the SPE was conditioned with methanol (MeOH) followed by water (H₂O). Subsequently, 400 µL of sample was mixed with 4 µL of 13C18-zearalenone (100 ng/mL), diluted with 396 µL of phosphate-buffered saline (PBS) before loading on the SPE plates. Plates were washed with H₂O and eluted with 2 x 200 µL of MeOH. 4 µL of the ISTD mix were added, and extracts were diluted with 396 µL of H₂O resulting in a total volume of 800 µL. Method details are provided in the SI.

UHPLC-MS/MS. LC separation was performed using a generic LC method (Jamnik et al., 2022; Preindl et al., 2019) with the detailed parameters presented in the SI and Table S4. In brief, an Agilent 1290 Infinity II UPLC system was equipped with an Acquity HSS T3 reversed-phase column (1.8 µm, 2.1 x 100 mm) and a VanGuard pre-column (1.8 µm), both purchased from Waters (Vienna). Retention times (RTs) for single compounds and potential RT shifts were determined by injecting standard solutions of single compounds or compound mixtures using concentrations of 0.5 ng/mL, 1 ng/mL, and 3 ng/mL.

Method optimization for multiple reaction monitoring (MRM) transitions included precursor ion (Q1) and product ion (Q3) selection, optimization of entrance potential (EP), collision energy (CE), and collision cell exit potential (CXP). Values of Q1, Q3, EP and CE were optimized for many compounds on a new Qtrap 7500 instrument (Sciex, Vienna) using method parameters of our previously published methods on a Qtrap 6500+ (Sciex, Vienna) as the basis (Braun et al., 2022; Hossain et al., 2024; Jamnik et al., 2022; Preindl et al., 2019). MRM parameters for all other analytes were

optimized using direct infusion of standard solutions (5–10 ng/mL) with a syringe pump at a flow rate of 7 μ L/min. Analytes were introduced either as mixtures or as individual solutions (see Table S5). Poorly ionizing compounds were re-analyzed individually and optimization was repeated. Two fragment ions were selected per compound as quantifier and qualifier ions. To ensure optimum data quality, dwell times were optimized for all targets to ensure enough datapoints per chromatographic peak. Fast polarity switching was used in scheduled MRM mode. More details, including ion source parameters and MRM transitions are listed in Table S6.

Quality control and quantification. A system suitability test (SST) was performed before and after each analytical sequence to track instrument performance (Table S7, and S8 and S9). External calibration with a multianalyte standard was used to quantify analytes and standards were prepared in 50% MeOH in H₂O (v/v, 1/1) at nine concentration levels (Table S10). Linear range and coefficient of determination (R^2) are reported in Table S11. Water was used as sample for the preparation of process blanks and 50% MeOH was used as solvent blank to control background concentrations and carryover. Three types of quality control (QC) samples were used including pooled samples (“non-spiked QC”), pooled samples spiked after extraction (“post-spiked QC”) and before extraction (“pre-spiked QC”). Furthermore, eleven isotopically labeled internal standards were used to track method performance during the analysis of 200 urine samples (details presented in the SI). Sciex OS software (version 4.1, Sciex, Vienna) was used to integrate peaks using the MQ4 algorithm, to build calibration curves and for quantification (more details presented in the SI). Further data evaluation and visualization was done in Microsoft Excel (v16.0; Microsoft Office), PowerPoint (v16.0; Microsoft Office), Origin 2021b (v 9.8; OriginLab Corporation), and Inkscape (version 1.2.2; Inkscape).

In-house validation. Following the European Commission Decision (EC) No. 2021/808 (Commission, 2021) performance parameters were evaluated over three days, including linearity, selectivity, matrix effects (signal suppression/enhancement, SSE), sensitivity (LOD and LOQ), trueness (extraction recovery, RE), intermediate precision (inter-day RSD, RSDR), and repeatability (intra-day RSD, RSDr). Sample preparation, measurements and data analysis were performed for each batch individually. Method linearity was evaluated by standards in solvent, and selectivity was tested by inspection of pooled urine and plasma samples for interferences. Matrix-matched calibration standards were prepared by spiking urine and plasma extracts at seven levels (Table S12). SSE (%) was calculated as the ratio of slopes of the matrix-matched and solvent calibration and was calculated for three batches (Matuszewski et al., 2003). LODs and LOQs were determined by analyzing spiked urine and plasma samples at low concentrations as replicates (n=9; see Table S13). LOD and LOQ were calculated as 3 and 10 times the standard deviation of the measured concentrations divided by the square root of the number of replicates (Cantwell, 2025).

RE was calculated based on determined concentrations in pre-spiked samples and using matrix matched calibration. RE was determined for two concentration levels, approximately 3 \times and 30 \times LOQ (n=3 each, see Table S14). Blank correction was performed and concentrations determined in non-spiked samples (n=4) was subtracted for RE calculation. RSDR and RSDr were determined as coefficient of variation (CV) of the RE within a single batch (n = 9) and within three validation batches. Full details of the validation approach are presented in the SI.

Validation criteria for trueness, precision, and repeatability based on EC guidelines. The EC criteria specifies performance criteria for analytical method validation. For trueness (RE), acceptable ranges depend on the spiked concentrations. For concentrations >10 μ g/kg, RE should be between 80% and 120%, 70%-120% is acceptable for concentrations between 1 μ g/kg and 10 μ g/kg,

whereas acceptable REs are in the range of 50%-120% for concentrations <1 μ g/kg. In analogy, acceptable values for RSDR and RSDr, determined as CV of RE, are also concentration-dependent. Following the guideline, CVs <16% are required for concentrations >1000 μ g/kg, while RSDs <22% is required for concentrations between 120 μ g/kg and 1000 μ g/kg. For concentrations <120 μ g/kg, the guideline recommends minimizing CVs for both RSDr and RSDR to technically feasible level. In the subsequent description, units of μ g/kg have been approximated to ng/mL, using the density of human urine/plasma as a basis for the conversion (more details presented in Table S15).

Conceptual framework of validation criteria in exposome-scale studies. While current EC criteria provide criteria for trueness, precision, and repeatability for analyte concentrations >120 ng/mL, a comprehensive framework for trace analysis in exposomic-scale investigations is still missing. However, the high number of analytes per assay and the required ultimate method sensitivity in the pg-ng/mL range warrant tailored adjustment of analytical figures of merit. To address this challenge, we propose novel tailored validation criteria based on empirical data. Validation criteria for exposomics were developed using published validation data of exposomics assays for large analyte panels (i.e., >90 analytes) covering a wide range of physicochemical properties (i.e., covering a logP range of >5) and results from our experiments as a rational basis (see Table S16) (Huang et al., 2023; Jamnik et al., 2022). Subsequently, the expected ranges for key parameters including RE, RSDR and RSDr values were derived as 5th and 95th percentiles of reported data and were established as new criteria for assessment and comparison in exposomic-type LC-MS investigations.

RESULTS AND DISCUSSION

Analyte selection. To ensure broad analyte coverage and to push the method to the edge of technical feasibility, diverse classes of environmental and food-related contaminants, other xenobiotics, and biotransformation products thereof were targeted. The analyte panel was based on our previous work and included PFAS chemicals, medical drugs, personal care products, air pollutants, pesticides, flame retardants, mycotoxins, industrial products, food processing by-products, plastic related chemicals, disinfectants and by-products, phytotoxins and phytoestrogens, and endogenous estrogens (Table S17). To complement the analyte panel, listed EDCs were selected to complement the analyte panel (Table S1). This includes the four disinfectants 2,4,5-trichlorophenol, 2,5-dichlorophenol, 2- and 4-phenylphenol, which are listed in the EPA's Endocrine Disruptor Screening Program (EPA|EDSP, www.epa.gov/comptox-tools). Previous studies reported low levels of these analytes in urine samples of adolescents and children with median concentrations of 0.02–4.64 ng/mL (Frederiksen et al., 2020; Guo et al., 2019). Ten additional personal care product-related compounds were selected based on EPA|EDSP, EPA's Ecotoxicology knowledgebase (EPA|ECOTOX) and the Comparative Toxicogenomics Database (CTD). Furthermore, 23 plastic related compounds, five bisphenols and 18 phthalates, were selected based on the EPA|EDSP, EPA|ECOTOX, Toxicity Values database (EPA|ToxValDB), EPA's Multimedia Monitoring Database (EPA|MMDB), and Chemicals of Emerging Concern (CEC) listed in the HBM4EU database (HBM4EU|CECscreen). These chemicals were previously reported in urine and serum from pregnant women, children and adolescents with median concentrations in the range of 0.09 – 46 ng/mL (Caballero-Casero et al., 2021; Frederiksen et al., 2020; Frigerio et al., 2020; Guo et al., 2021; Gys et al., 2021; Jala et al., 2022; Varghese et al., 2022). Five flame retardants, tri-n-butyl phosphate, triphenyl phosphate, tris(2-chloroethyl)phosphate, bis(1,3-dichloro-2-propyl)phosphate, and bis(1-chloro-2-propyl)phosphate, were selected based on EPA|EDSP and EPA|ECOTOX lists and due to reports in human urine and blood samples (Gao et al., 2020; Li et al., 2020; Siddique et al., 2020).

Finally, five PFAS substances (perfluorobutane sulfonic acid, perfluorodecanoic acid, perfluorohexanoic acid, perfluorononanoic acid, and perfluoroundecanoic acid), three drugs (duloxetine, tramadol and venlafaxine), three bioactive phytochemicals (chrysin, chalcone, and kaempferol), and five pesticides (glyphosate, vinclozolin, propiconazole, tebuconazole, and metribuzin) listed in EPA ECOTOX or EPAEDSP database were added.

LC-MS method optimization. To ensure favorable method performance, RT windows and dwell times for all scheduled MRM transitions were optimized based on observed peak widths and RT shifts in biological matrices. Optimization of window width and dwell times is essential to balance total cycle time, number of datapoints per peak, RT stability and the method sensitivity. Based on empirical observation of the chromatographic peak widths, method development aimed for a total cycle time <850 ms which could be achieved for the majority of all analytes (see Figure 1b). Due to unstable RTs in matrix samples, diethyl dithiophosphate, ibuprofen, perfluorodecanoic acid and perfluorononanoic acid, as well as six phthalates were monitored during the entire runtime. After method optimization, the number of datapoints per peak was investigated in spiked samples ($n=3$) to ensure sufficient data quality for quantification with 91% of the analytes having >7 datapoints per peak (Table S18 and Figure S2(B)) which can be regarded as sufficient for quantification (Zeng & Bateman, 2023).




Figure 1. Scheduled Multiple Reaction Monitoring (MRM) windows (A) and MRM cycle time (ms, B) as function of the chromatographic retention time (min) of the optimized LC-MS/MS method. Extraction recovery (RE, %, C) and signal suppression and enhancement (SSE, %, D) in dependence of chromatographic retention time (min) for analytes in urine and plasma. Four analytes (2,5-dichlorophenol, xanthohumol, bisphenol M, and triclosan) are not shown in panel (C) due to high RE caused by interferences and/or matrix effects.

Development of tailored validation criteria for exposomics. Several different guidelines for analytical method validation provide valuable frameworks of criteria for method validation in diverse applications. The choice of the used criteria typically depends on the specific research area. Among them, the EC guidelines provide a comprehensive set of criteria (EC criteria) for various application fields, however, they do not address key challenges inherent to exposome-scale LC-MS/MS analysis. In particular, the methodological complexity and the extreme sensitivity requirements needed to detect large mixtures of highly diverse analytes at trace levels within a single assay is not covered. To address this gap, tailored validation criteria were developed and are proposed herein, based on published validation data from a set of exposomics-scale LC-MS/MS methods quantitatively assessing >90 analytes and

spanning a broad logP range of >5. Together with the data presented herein, these results were used as the basis for the implementation of new criteria for method validation of exposomics-scale multianalyte assays using 5th and 95th percentiles of empirical validation data. Figure 2A-C show the distribution of RE and RSD reported in two published methods (Huang et al., 2023; Jamnik et al., 2022) and results from the study at hand. The 5th and 95th percentiles of RE were calculated as 42% and 134%, respectively, highlighting certain limitations in terms of analyte recovery for currently used LC-MS workflows in case of sets of extremely diverse analytes at trace concentrations. These empirically determined limits were subsequently proposed as the acceptable range for RE in quantitative exposomics (Figure 2A). In analogy, acceptable limits for assay repeatability and precision were determined as $RSDr < 37\%$ and $RSDr < 42\%$ (Figure 2B and Figure 2C). These custom validation criteria were then used to complement the EC criteria for the relevant concentration range in exposome research (i.e., <1 ng/mL), and were integrated into our proposed method validation framework (see Figure 2D and Table S15). Figure 2 illustrates a schematic representation of the criteria provided in the EC guideline and proposed new data-derived criteria for exposomics including relevant concentration ranges. Furthermore, a simple three-level scheme is proposed to categorize targeted analytes based on validation results: quantitative data is reported as data quality category A, semi-quantitative data is reported as data quality category B, and qualitative data as data quality category C. Detailed definitions for A, B and C are described below (section *Classification of data quality*)

Results of the in-house validation. Excellent method linearity was achieved, with 88% of compounds achieving R^2 values exceeding 0.99 in pure solvents (see Table S11). Sensitivity was also excellent for the majority of analytes with LOD values below 1.0 ng/mL observed for 73% of compounds in urine and 74% in plasma. Approximately 10% of the analytes showed exceptionally low LODs (< 0.001 ng/mL) in both sample matrices (Figure 2E). For approximately 70% of the analytes, RE values were in the range of 42%–134% in both sample matrices (Figure 3A) and 80% and 77% of the analytes showed $RSDr < 37\%$ in plasma and urine, respectively. Notably, 63% displayed $RSDr$ values below 20% (Figure 3B). $RSDr$ was determined with 67% of all compounds showing $RSDr < 37\%$, and 38% and 57% of compounds performing $RSDr < 20\%$ in plasma and urine, respectively (Figure 3C). Matrix effects were evaluated and effective sample cleanup with SPE was demonstrated as previously reported for a smaller subset of the analyte panel (Gu et al., 2025). In total, 61% of the analytes in urine and 70% in plasma showed SSE within 60%–140% (Figure 3D). Full validation results are provided in Table S14 and Table S19.

Classification of data quality: requirements for quantification, semi-quantification and qualification in quantitative exposomics. To be in line with the EC criteria, analytes need to meet the following requirements for at least one spiking level to be reported as quantitative data (category A): 1.) RE within 50%–120% for concentrations ≤ 1 ng/mL, between 70%–120% for concentrations between 1–10 ng/mL and 80%–120% for concentrations between 10–120 ng/mL or higher. 2.) $RSDr$ and $RSDr \leq 22\%$ are required for concentrations within 120–1000 ng/mL, $RSDr$ and $RSDr \leq 16\%$ for concentrations >1000 ng/mL, whereas no clear limit is defined for $RSDr$ and $RSDr$ for compounds at trace and ultra-trace levels (see Figure 2D). Using the tailored criteria, data was reported as quantitative if the following requirements were met for at least one spiking level: 1.) RE in the range of 42%–134% and 2.) $RSDr \leq 37\%$ and $RSDr \leq 42\%$, respectively (see Figure 2D). To implement a clear scheme for reporting quantitative, semi-quantitative and qualitative data in exposomics, two additional categories were introduced for the classification of the data quality: data was classified as semi-quantitative (category B) if either RE or RSD values failed to meet quality criteria for quantitative data described above. Qualitative data (category C) was reported for compounds




Figure 2. Distribution of extraction recovery (RE, %, A), repeatability (RSDr, %, B), and precision (RSDR, %, C) reported in two previously published exposome-scale targeted LC-MS/MS studies (Huang et al., 2023; Jamnik et al., 2022) and data generated in the study at hand. Empirical data was used to derive tailored validation criteria for exposome-type method validation to complement and extend EC criteria (EC No. 2021/808) for extremely low-abundance analyte concentrations in multianalyte assays (D). Panel (E) shows the distribution of the LOD values in the present study, highlighting the need to extend the considered concentration range of validation guidelines to <1 ng/mL when measuring toxic and bioactive chemicals. NA, not available data, i.e., LODs not determined.

with stable RT and signal intensity, but accuracy and repeatability could not be evaluated due to poor sensitivity or high background.

In Figure 3F the number of analytes meeting data quality criteria in the EC guidelines and using our newly proposed approach are compared. Using the EC criteria, 170 analytes were classified as “quantitative” in urine, whereas application of the tailored criteria results in 163 quantifiable analytes. In plasma, the EC criteria resulted in 162 quantifiable analytes, compared to 159 using our tailored criteria. Noteworthy, most analytes (155 in urine and 146 in plasma) met the quantification requirements of both validation schemes. A subset of compounds (15 in urine and 16 in plasma) met quantification level A under the EC criteria, yet were classified as level B when using the tailored criteria, primarily due to relatively high RSDr and RSDr values at concentrations below 120 ng/mL, which is effectively not defined in the EC criteria. This difference reflects the EC criteria's greater tolerance for poor precision and repeatability at extremely low analyte concentrations. In contrast, another subset of analytes (8 in urine and 13 in plasma) was classified as semiquantitative (category B) using EC criteria but as quantitative (category A) using the tailored classification scheme due to a higher tolerance in regards of RE. Here, the tailored criteria allow higher tolerance for trueness, due to the complexity of large-scale exposome-scale assay. This takes into consideration that correction for “trueness” (i.e., lower RE) is often possible if the method's precision is acceptable and a suitable correction approach (e.g., using internal standards

or quality control samples) is used. Given the complementary nature of both approaches and the complex nature of large-scale multi-analyte LC-MS/MS assays, compounds meeting quality criteria of either of the validation schemes were considered as suitable for quantification. In total, 76% and 75% of all analytes met requirements for quantification (category A) in urine and plasma, respectively. Another 10% and 15% of all analytes met requirements for semi-quantification (category B), whereas the remaining 14% and 10% can only be reported as qualitative data (i.e., detected and non-detected, category C) in urine and plasma, respectively (see Figure 3E).

Benchmarking coverage, sensitivity, and throughput. The presented and fully validated exposomics workflow allows analyzing diverse exposure chemicals covering a wide range of polarities and concentrations. To benchmark its performance, method performance and coverage were compared to three relevant published LC-MS/MS protocols (see Table S16 for full details). Noteworthy, we could present the dataset covering the widest range of analyte polarities and achieved comparable or better sensitivity for the majority of analytes (Figure 4). Especially, the fraction of extremely sensitive analytes (LOQ <0.0033 ng/mL) was approximately two times higher than in previously published studies (see Figure 4B) (Huang et al., 2023; Jamnik et al., 2022; Zhu et al., 2021).

While the absolute number of analytes was highest in the method presented by Huang et al. (2023) (n=325), the presented

workflow includes the most diverse set of analytes (n=234) in terms of analyte polarity. This is a critical metric to show sufficient coverage and scalability potential for exposomics applications and highlights competitive coverage with other state-of-the-art LC-




Figure 3. Results of the in-house method validation based on performance: analytes with RE (%), RSDr (B, %), RSDR (C, %), and SSE (%), within the range required for quantitative exposomics. Classified validation results following EC criteria and proposed tailored validation criteria using three categories: quantification, semi-quantification, and qualification (E). Number of analytes meeting quantification criteria for at least one spiking level following the European Commission guideline (EC, No. 2021/808) and the proposed tailored set of criteria (F). The results show that the majority of quantitative analytes in the present study fully comply with the established EC criteria.

MS/MS-based workflows for exposomics (Figure 4A) (Huang et al., 2023; Jamnik et al., 2022; Zhu et al., 2021). Clearly, it demonstrates further possibilities to scale up the developed workflow. Assay sensitivity was benchmarked by comparing method LOQs with other relevant literature and a special focus on analytes with good method sensitivity and LOQs <3.3 ng/mL. In general, this number was comparable to other assays, however, outstanding method sensitivity was demonstrated for 10% of the assay compounds with LOQs <0.0033 ng/mL in both, urine and plasma. In contrast, only 4%-7% of the analytes had comparable LOQs in other published work (see Figure 4B and Table S16). The applicability of the developed method for high-throughput exposomics was assessed and compared to other protocols in terms of sample requirements, solvent consumption and total analysis time (Table S20). Using the SPE method presented herein, 1.4 mL of solvents are required per sample, compared to 0.8 mL estimated for a previously published PPT protocol (Jamnik et al., 2022). Similar calculations were performed for other published work resulting in even higher solvent consumption of approximately 6 mL per sample calculated for the sample preparation protocol of Huang et al. (2023). While these estimations include some uncertainties, this further highlights the high efficiency of our workflow which would result in a theoretical reduction of 4.6 L solvents calculated for the extraction of 1000 samples. Finally,

sample preparation in 96-well plates is time efficient. It was estimated that the total required laboratory time of 3 days for processing 1000 samples including SPE cleanup is sufficient whereas the same number of samples would need roughly 30 working days using a traditional PPT protocols in Eppendorf tubes (Gu et al., 2025; Jamnik et al., 2022).




Figure 4. (A) Coverage of the chemical space of multi-analyte exposomics assays in terms of octanol-water partition coefficient (LogP) values for analytes in this study (n=234) and three published papers: Zhu et al. (2021; n=121), Jamnik et al. (2022; n=96), and Huang et al. (2023; n=325). (B) Limit of quantification (LOQ, ng/mL) for analytes in plasma (n=206) and urine (n=208) in this study, in urine (n=114) in Zhu et al. (2021), serum (n=92) and urine (n=87) in Jamnik et al. (2022), and plasma (n=295) in Huang et al. (2023). Note: Due to the lack of reported LOQ values in Zhu et al. (2021), we estimated LOQs as 3.33x the reported LODs in (B); The LOQ value of di-2-propylheptyl phthalate (DHPH) in this study was excluded from panel B due to high background levels in the non-spiked samples. LogP values listed in PubChem were used for panel (A).

Real-life exposure patterns in pregnant women from the YPOPS cohort. To demonstrate favorable sensitivity and general performance of the new workflow and its applicability to real-life research questions, 200 urine samples of women (n=50) from a US pregnancy cohort were analyzed. Pregnant women present a highly vulnerable population group and expected exposure levels are typically lower than in the average population. Yet, 143 out of 234 target analytes were detected with 37 exposure compounds present with very high detection frequencies (DF) above 70% (Figure 5). The detected and quantified biomarkers of exposure included chemicals from various exposure routes, application fields and compound classes including antibiotics, mycotoxins, medical drugs, pesticides, personal care products, phytotoxins, phytoestrogens, plasticizers, and PFAS chemicals.

Examples of selected toxicants are briefly discussed below. The very data-rich detailed results including the classification of the data quality based on the obtained validation data is presented in Tables S21 and S22. A comprehensive discussion of the comprehensive exposure patterns is beyond the scope of this work.

Urinary cotinine, a biomarker for cigarette smoke, was determined in 60% of the samples (median 0.22 ng/mL, range 0.1-23 ng/mL), demonstrating again the high sensitivity with a LOD

value of 0.03 ng/mL. A phase-I-metabolite of cotinine, trans-3-hydroxycotinine, was quantified in 74% of the samples at a lower median concentration (median 0.16 ng/mL, range 0.07-25 ng/mL). Importantly, these levels typically indicate environmental smoke exposure rather than active smoking.

Several PFAS chemicals including perfluorooctanoic acid (PFOA; median 0.024, range 0.024 – 0.12) and perfluorooctanesulfonic acid (PFOS; median 0.0065, range 0.0035 – 0.014) were quantified in up to 67% of the samples. Typically, PFASs are assessed in serum for HBM purpose, however, highly sensitive approaches are also able to determine them in urine samples (Zhang et al., 2015).

Nine out of 12 monitored bisphenol derivatives or metabolites were quantified. As expected, BPA-glucuronide was by far the most frequently detected analyte of this class, however, only level C data quality (qualitative results) was generated. Accordingly, the parent molecule BPA was frequently detected (73%) (median 0.11 ng/mL, range 0.09 – 0.53 ng/mL). For bisphenols other than BPA, no glucuronide or sulfate conjugates were directly accessed, consequently, the comparison of total exposure levels is not straight forward. However, these results clearly demonstrate the feasibility of this broad method to also accurately and directly quantify highly polar xenobiotic metabolites such as glucuronides.

The pesticides, acetamiprid, imidacloprid (Figure 5), and the pesticide metabolite 2-isopropyl-4-methyl-6-hydroxypyrimidine (IMPY) were detected in most samples. The results were compared to other published datasets of pesticide exposure in pregnant women from Asia and Europe. For acetamiprid and imidacloprid, DFs of 91% and 96% were observed in our study with concentration levels ranging <0.00023-0.044 ng/mL and 0.025-1.8 ng/mL, respectively. Lower DFs (18% and 26%) and comparable urinary concentrations (<0.0025-1.3 ng/mL and <0.025-0.95 ng/mL) were reported for 617 urine samples collected from 62 pregnant women in Japan (Suwannarin et al., 2024). Similarly, IMPY, which is a key metabolite of the organophosphate pesticide diazinon, was detected at concentration levels between <0.0086-2.1 ng/mL (DF=72%) in the present study. In a study from Spain, IMPY was present at concentrations between <1.6-744.2 ng/mL and was detected in 12% of urine samples from pregnant women (n=573) in Valencia (Llop et al., 2017).

The data obtained from this high-performance multi-analyte assay also demonstrated proper performance when comparing it to a more tailored HBM assay covering multiple mycotoxins that was applied to the same set of sample before (Krausová et al., 2024). The most concentrated and prevalent fungal toxin, deoxynivalenol was determined in 50% of the samples at a median concentration of 3.2 ng/mL (Table S21), notably in its native form. The median in the tailored HBM assay that measured total

deoxynivalenol after enzymatic deconjugation of glucuronide and sulfate conjugates, was reported as 23 ng/mL. Given the fact that the vast majority of deoxynivalenol is present conjugated in human urine (Turner et al., 2011; Warth et al., 2013), this demonstrated comparable performance. The same is true for the *Alternaria* toxin alternariol. Alternariol-monomethylether were in line with quantitative results from this cohort analyzed before and a cohort of primary school children in Austria (Ayeni et al., 2023). Ochratoxin A demonstrated very similar results with reported median concentrations of 0.21 ng/mL (this study) and 0.2 ng/mL (Krausová et al., 2024). For citrinin, a nephrotoxic mycotoxin, the sensitivity of the presented assay was even enhanced, so that a low number of samples indicated exposure while the previously published, tailored HBM assay was unable to detect this food contaminant in urine. The same was true of enniatin B, a *Fusarium* toxin, that was quantified in a low number of samples at 0.0019 ng/mL. Beauvericin was quantified in human urine for the first time to the best of our knowledge (Table S21).

In addition, many other natural toxins and bioactives have been identified and quantified. Some of them are reported in the urine of US citizens for the first time to the best of our knowledge, while others, especially those occurring at higher concentrations, have been reported before (Warth et al., 2017). The first time identifications include Aristolactam I, a nephrotoxic plant toxin and scopolamine, a tropane alkaloid that was recently quantified in pooled US serum (NIST SRM 1958; Verri Hernandez & Warth, 2024) but not yet in urine.

The exposome is defined as the totality of the exposure to the (chemical) environment over the lifespan and longitudinal studies are of special relevance. In the investigated sample cohort of the YPOPS study, 50 pregnant women were sampled at four timepoints during pregnancy. Figure 6 illustrates observed time trends of urinary concentrations for ten biomarkers across the four gestational periods. The trends included a continuous increase of a monitored hormone (estradiol), an increase across the first three periods followed by a decrease in the final period (hormone metabolite, estradiol-3-sulfate), stable (pesticide, 2,4 dichlorophenoxyacetic acid), and decrease across the first three periods followed by an increase in the final period for diet related product (enterolactone), drug metabolite (ibuprofen glucuronide), personal care product (methylparaben), perfluorooctanoic acid, plastic related product (mono-(2-ethyl-5-hydroxyhexyl) phthalate and bisphenol A glucuronide), and UV filter (benzophenone 1). More detailed illustrations for the urinary concentrations are presented in Figure 5. Detailed results and descriptive statistics for all analytes are summarized in Table S22.




Figure 5. Concentrations of frequently detected compounds (DF>70% and data quality level A) in urine samples from pregnant women in the YPOPS cohort. Note: Concentrations below the limit of quantification (LOQ) were imputed with LOQ/2. 2,4-D, 2,4 Dichlorophenoxyacetic acid; GLCA, glucuronide; 5-OH-MEHP, Mono-(2-ethyl-5-hydroxyhexyl) phthalate; MCMHP, Mono-[(2-carboxymethyl) hexyl] phthalate; MECPP, Mono-(2-ethyl-5-carboxypentyl) phthalate; MEOHP, Mono-(2-ethyl-5-oxohexyl) phthalate.

1




Figure 6. Longitudinal concentration profiles of frequently detected compounds (DF>70%) in urine samples from 50 pregnant females across four gestational periods in the Yale Pregnancy Outcome Prediction Study (YPOPS). Some exposures are highly stable throughout pregnancy while others are dynamic and exhibit changing concentration levels. Here, concentrations below the limit of quantification (LOQ) were replaced with a value of LOQ/2. PFAS, per- and polyfluoroalkyl substances; 2,4-D, 2,4 dichlorophenoxyacetic acid; GLCA, glucuronide; 5-OH-MEHP, mono-(2-ethyl-5-hydroxyhexyl) phthalate.

LIMITATIONS

While the robustness, sensitivity as well as time and cost-efficiency of the developed workflow were clearly demonstrated, several limitations remain. Although the presented analyte panel was well designed and aims for high diversity and representativeness, still only a limited subset of the chemical exposome can be covered in a single assay. For a comprehensive assessment of the human exposome, even wider chemical coverage or combinations of complementary methods are required. Furthermore, several very polar compounds have been lost during the SPE process due to weak retention with the sorbent and substantial matrix effects were observed for several analytes. Combinations of different orthogonal approaches might address these limitations, yet it would be time-consuming and costly. Clearly, broad coverage exposome-type assays need to find a well-balanced compromise between analyte coverage, efficient sample cleanup, assay sensitivity and time requirements. The overall performance of the method in the real-life case study yielded highly information-rich exposure patterns, indicating that the chosen pragmatic approach seems fit-for-purpose.

Traditionally used criteria for analytical method validation show limited applicability for exposome research because of the required ultimate sensitivity and the method complexity. Our newly established, validation criteria aim to propose a conceptually novel and evidence-based framework to systematically extend established guidelines for analytical method validation to meet requirements of quantitative exposomics. This is of great relevance to ensure comparable reporting standards and method performance. Yet, our proposed tailored validation criteria rely on data derived from only a small number of published exposome-scale studies (n=3) due to limited availability of exposome-

scale LC-MS/MS assays. Further integration of additional datasets would be needed before claiming to provide a comprehensive framework for method validation in exposomics-scale HBM.

CONCLUSION AND OUTLOOK

In summary, a high-performance workflow for exposome analysis characterized by broad chemical coverage and excellent method sensitivity across diverse human matrices was developed. The analytical pipeline was validated based on existing guidelines and using a conceptually new and tailored framework of validation criteria designed for exposomics. This is necessary to face typically observed challenges in exposomics, such as broad exposome-scale analyte coverage and compounds present at trace concentrations. Method performance was characterized and offers a robust foundation for quantifying numerous priority exposures in population-scale studies. The versatile assay extends analyte coverage beyond traditional food and environmental toxicants and includes biotransformation products, medical drugs and microbiome-related compounds besides typically monitored toxins and CECs. Sample throughput presents one of the critical bottlenecks in large-scale exposomics and we could reduce the required laboratory work by approximately a factor of ten using SPE in 96-well plates. The scalability of the method was demonstrated and now offers the required capacities to perform large-scale population-based studies for ExWAS. Such large-scale studies are finally needed to enhance our understanding of the chemical exposome, which is prerequisite to improve disease prevention in public health and personalized medicine. Due to the excellent method sensitivity and sample throughput, our validated workflow offers tools to increase our understanding of the chemical exposome and will be used for in large-scale population studies. Such studies will be crucial for unraveling the

relationships between (co-)exposure to diverse environmental chemicals and human health in the future.

ASSOCIATE CONTENT

Supplementary information (SI) provided in the PDF Document is included and provided with additional information regarding the method and the results.

ETHICS DECLARATIONS

Competing interests

The authors declare no competing interests.

ACKNOWLEDGMENT

The authors would like to thank all members of Warth Laboratory for providing valuable feedback and discussions, especially Maren Kirchner for her assistance during sample pretreatment. They express their gratitude to the Mass Spectrometry Center of the Faculty of Chemistry at the University of Vienna for technical support during the measurements. Moreover, the authors greatly acknowledge the Yale University Reproductive Sciences (YURS) Biobank (HIC#1309012696), especially Jane O'Bryan, Lauren Perley, and Michelle Silasi and all the women who donated samples.

This work was supported by the University of Vienna, the China Scholarship Council (CSC), the Austrian Science Fund (FWF, P33188), the European Research Council, ERC, (EXPOMET, 101043321), and Exposome Austria, the national Node of the EIRENE research infrastructure. Views and opinions expressed are those of the authors only and do not necessarily reflect those of the ERC or CSC Executive Agency. None of the granting authorities can be held responsible for them.

REFERENCES

- Ayeni, K. I., Jamnik, T., Fareed, Y., Flasch, M., Braun, D., Uhl, M., Hartmann, C., & Warth, B. (2023). The Austrian children's bio-monitoring survey 2020 part B: Mycotoxins, phytotoxins, phytoestrogens and food processing contaminants. *Food and Chemical Toxicology*, 182, 114173.
- Braun, D., Abia, W. A., Šarkanj, B., Sulyok, M., Waldhoer, T., Erber, A. C., Krska, R., Turner, P. C., Marko, D., Ezekiel, C. N., & Warth, B. (2022). Mycotoxin-mixture assessment in mother-infant pairs in Nigeria: From mothers' meal to infants' urine. *Chemosphere*, 287(Pt 2), 132226. <https://doi.org/10.1016/j.chemosphere.2021.132226>
- Braun, D., Schernhammer, E., Marko, D., & Warth, B. (2020). Longitudinal assessment of mycotoxin co-exposures in exclusively breastfed infants. *Environment International*, 142, 105845.
- Caballero-Casero, N., Castro, G., Bastiaensen, M., Gys, C., van Larebeke, N., Schoeters, G., & Covaci, A. (2021). Identification of chemicals of emerging concern in urine of Flemish adolescents using a new suspect screening workflow for LC-QTOF-MS. *Chemosphere*, 280, 130683. <https://doi.org/10.1016/j.chemosphere.2021.130683>
- Cai, J., & Yan, Z. (2021). Re-Examining the Impact of Minimal Scans in Liquid Chromatography–Mass Spectrometry Analysis. *Journal of the American Society for Mass Spectrometry*, 32(8), 2110–2122.
- Cantwell (ed.) Eurachem Guide: The Fitness for Purpose of Analytical Methods: A Laboratory Guide to Method Validation and Related Topics: Third edition (2025), (2025). <http://www.eurachem.org>
- Commission, E. (2021). Commission Implementing Regulation (EU) 2021/808 of 22 March 2021 on the performance of analytical methods for residues of pharmacologically active substances used in food-producing animals and on the interpretation of results as well as on the methods to be used for sampling and repealing Decisions 2002/657/EC and 98/179/EC. *Off. J. Eur. Union*, 180, 84–109.
- Frederiksen, H., Nielsen, O., Koch, H. M., Skakkebaek, N. E., Juul, A., Jørgensen, N., & Andersson, A.-M. (2020). Changes in urinary excretion of phthalates, phthalate substitutes, bisphenols and other polychlorinated and phenolic substances in young Danish men; 2009–2017. *International journal of hygiene and environmental health*, 223(1), 93–105.
- Frigerio, G., Campo, L., Mercadante, R., Santos, P. M., Missineo, P., Polledri, E., & Fustinoni, S. (2020). Development and validation of a liquid chromatography/tandem mass spectrometry method to quantify metabolites of phthalates, including di-2-ethylhexyl terephthalate (DEHTP) and bisphenol A, in human urine. *Rapid Communications in Mass Spectrometry*, 34(13), e8796.
- Gao, D., Yang, J., Bekele, T. G., Zhao, S., Zhao, H., Li, J., Wang, M., & Zhao, H. (2020). Organophosphate esters in human serum in Bohai Bay, North China. *Environmental Science and Pollution Research*, 27, 2721–2729.
- Gu, Y., Feuerstein, M. L., & Warth, B. (2025). High-Throughput Solid Phase Extraction for Targeted and Nontargeted Exposomics. *analytical Chemistry*. <https://doi.org/10.1021/acs.analchem.4c06177>
- Gu, Y., Peach, J. T., & Warth, B. (2023). Sample preparation strategies for mass spectrometry analysis in human exposome research: Current status and future perspectives. *TrAC Trends in Analytical Chemistry*, 117151.
- Guo, C., Zhao, X., Jin, J., Wang, L., Tan, D., Chen, J., & Ni, Y. (2021). The dose effect of dansyl chloride on the derivative products of bisphenols and its application for the determination of bisphenols in human serum by high-performance liquid chromatography–tandem mass spectrometry. *Journal of Separation Science*, 44(16), 3052–3060.
- Guo, J., Wu, C., Zhang, J., Jiang, S., Lv, S., Lu, D., Qi, X., Feng, C., Liang, W., & Chang, X. (2019). Anthropometric measures at age 3 years in associations with prenatal and postnatal exposures to chlorophenols. *Chemosphere*, 228, 204–211.
- Gys, C., Bastiaensen, M., Malarvannan, G., Bamai, Y. A., Araki, A., & Covaci, A. (2021). Short-term variability of bisphenols in spot, morning void and 24-hour urine samples. *Environmental Pollution*, 268, 115747.
- Hossain, M. Z., Feuerstein, M. L., Gu, Y., & Warth, B. (2024). Scaling up a targeted exposome LC-MS/MS biomonitoring method by incorporating veterinary drugs and pesticides. *Analytical and Bioanalytical Chemistry*, 1–14.
- Huang, W., Pan, X.-F., Tang, S., Sun, F., Wu, P., Yuan, J., Sun, W., Pan, A., & Chen, D. (2023). Target Exposome for Characterizing Early Gestational Exposure to Contaminants of Emerging Concern and Association with Gestational Diabetes Mellitus. *Environmental science & technology*, 57(36), 13408–13418.
- Jagani, R., Pulivarthi, D., Patel, D., Wright, R. J., Wright, R. O., Arora, M., Wolff, M. S., & Andra, S. S. (2022). Validated single urinary assay designed for exposomic multi-class biomarkers of common environmental exposures. *Analytical and bioanalytical chemistry*, 414(19), 5943–5966.
- Jala, A., Varghese, B., Dutta, R., Adela, R., & Borkar, R. M. (2022). Levels of parabens and bisphenols in personal care products and urinary concentrations in Indian young adult women: implications for human exposure and health risk assessment. *Chemosphere*, 297, 134028.
- Jamnik, T., Flasch, M., Braun, D., Fareed, Y., Wasinger, D., Seki, D., Berry, D., Berger, A., Wisgrill, L., & Warth, B. (2022). Next-generation biomonitoring of the early-life chemical exposome in neonatal and infant development. *Nat Commun*, 13(1), 2653. <https://doi.org/10.1038/s41467-022-30204-y>
- Koual, M., Tomkiewicz, C., Cano-Sancho, G., Antignac, J.-P., Bats, A.-S., & Coumoul, X. (2020). Environmental chemicals, breast

- cancer progression and drug resistance. *Environmental Health*, 19, 1-25.
- Krausová, M., Ayeni, K. I., Gu, Y., Borutzki, Y., O'Bryan, J., Perley, L., Silasi, M., Wisgrill, L., Johnson, C. H., & Warth, B. (2024). Longitudinal biomonitoring of mycotoxin exposure during pregnancy in the Yale pregnancy Outcome Prediction study. *Environment International*, 109081.
- Lai, Y., Koelmel, J. P., Walker, D. I., Price, E. J., Papazian, S., Manz, K. E., Castilla-Fernández, D., Bowden, J. A., Nikiforov, V., & David, A. (2024). High-resolution mass spectrometry for human exposomics: Expanding chemical space coverage. *Environmental Science & Technology*, 58(29), 12784-12822.
- Li, Y., Li, D., Chen, J., Zhang, S., Fu, Y., Wang, N., Liu, Y., & Zhang, B. (2020). Presence of organophosphate esters in plasma of patients with hypertension in Hubei Province, China. *Environmental Science and Pollution Research*, 27, 24059-24069.
- Lin, J., Ding, J., Di, X., Sun, W., Chen, H., & Zhang, H. (2022). Association between prenatal antibiotics exposure and measures of fetal growth: A repeated-measure study. *Ecotoxicology and Environmental Safety*, 244, 114041.
- Llop, S., Murcia, M., Iñiguez, C., Roca, M., González, L., Yusà, V., Rebagliato, M., & Ballester, F. (2017). Distributions and determinants of urinary biomarkers of organophosphate pesticide exposure in a prospective Spanish birth cohort study. *Environmental Health*, 16, 1-15.
- Manz, K. E., Feerick, A., Braun, J. M., Feng, Y.-L., Hall, A., Koelmel, J., Manzano, C., Newton, S. R., Pennell, K. D., & Place, B. J. (2023). Non-targeted analysis (NTA) and suspect screening analysis (SSA): a review of examining the chemical exposome. *Journal of Exposure Science & Environmental Epidemiology*, 33(4), 524-536.
- Matuszewski, B. K., Constanzer, M., & Chavez-Eng, C. (2003). Strategies for the assessment of matrix effect in quantitative bioanalytical methods based on HPLC- MS/MS. *analytical Chemistry*, 75(13), 3019-3030.
- Preindl, K., Braun, D., Aichinger, G., Sieri, S., Fang, M., Marko, D., & Warth, B. (2019). A Generic Liquid Chromatography-Tandem Mass Spectrometry Exposome Method for the Determination of Xenoestrogens in Biological Matrices. *Anal Chem*, 91(17), 11334-11342. <https://doi.org/10.1021/acs.analchem.9b02446>
- Siddique, S., Harris, S. A., Kosarac, I., Latifovic, L., & Kubwabo, C. (2020). Urinary metabolites of organophosphate esters in women and their relationship with serum lipids: An exploratory analysis. *Environmental Pollution*, 263, 114110.
- Suwannarin, N., Isobe, T., Nishihama, Y., Ito, Y., Kamijima, M., Ebara, T., Sugiura-Ogasawara, M., Nishikawa, N., Nakai, K., & Minamikawa, Y. (2024). Concentrations of neonicotinoid insecticides and their metabolites in multiple urine samples collected from pregnant women in Japan. *Environmental Research*, 240, 117506.
- Turner, P. C., Hopton, R. P., White, K. L., Fisher, J., Cade, J. E., & Wild, C. P. (2011). Assessment of deoxynivalenol metabolite profiles in UK adults. *Food and Chemical Toxicology*, 49(1), 132-135.
- Varghese, B., Jala, A., Das, P., Borkar, R. M., & Adela, R. (2022). Estimation of parabens and bisphenols in maternal products and urinary concentrations in Indian pregnant women: daily intake and health risk assessment. *Environmental Science and Pollution Research*, 1-14.
- Verri Hernandes, V., & Warth, B. (2024). Bridging targeted (Zeno MRM-HR) and untargeted (SWATH) LC-HRMS in a single run for sensitive exposomics. *Analytical Chemistry*, 96(31), 12710-12717.
- Warth, B., Spangler, S., Fang, M., Johnson, C. H., Forsberg, E. M., Granados, A., Martin, R. L., Domingo-Almenara, X., Huan, T., & Rinehart, D. (2017). Exposome-scale investigations guided by global metabolomics, pathway analysis, and cognitive computing. *Analytical Chemistry*, 89(21), 11505-11513.
- Warth, B., Sulyok, M., Berthiller, F., Schuhmacher, R., & Krska, R. (2013). New insights into the human metabolism of the Fusarium mycotoxins deoxynivalenol and zearalenone. *Toxicology letters*, 220(1), 88-94.
- Wild, C. P. (2005). Complementing the genome with an "exposome": the outstanding challenge of environmental exposure measurement in molecular epidemiology. *Cancer Epidemiol Biomarkers Prev*, 14(8), 1847-1850. <https://doi.org/10.1158/1055-9965.Epi-05-0456>
- Zeng, W., & Bateman, K. P. (2023). Quantitative LC-MS/MS. 1. Impact of Points across a Peak on the Accuracy and Precision of Peak Area Measurements. *Journal of the American Society for Mass Spectrometry*, 34(6), 1136-1144.
- Zhang, T., Sun, H., Qin, X., Gan, Z., & Kannan, K. (2015). PFOS and PFOA in paired urine and blood from general adults and pregnant women: assessment of urinary elimination. *Environmental Science and Pollution Research*, 22, 5572-5579.
- Zhu, H., Chinthakindi, S., & Kannan, K. (2021). A method for the analysis of 121 multi-class environmental chemicals in urine by high-performance liquid chromatography-tandem mass spectrometry. *J Chromatogr A*, 1646, 462146. <https://doi.org/10.1016/j.chroma.2021.462146>